


3-level balloon dilatator

For single use

3-level balloon dilatator for single use. The balloon diameter can be set to 3 different levels. The balloon length is 50 mm and the catheter length is 230 cm.



- For insertion in **all working channels of diameter 2.8 mm and larger**
- Balloon diameter **can be set to 3 levels**
- **Individual, sterile** packaging
- **Two techniques with a single balloon:**
 - As TTS balloon with the help of an embedded mandrin
 - As OTW balloon via guidewire .035" after removing the mandrin

➤ *You can find our standard sizes on the back of this page.*

Esophageal Balloon Dilator

Part number	Catheter length	Catheter Ø	Balloon Ø			Balloon length
34106PRO	230 cm	7 Fr	6 mm	7 mm	8 mm	55 mm
34108PRO	230 cm	7 Fr	8 mm	9 mm	10 mm	55 mm
34110PRO	230 cm	7 Fr	10 mm	11 mm	12 mm	55 mm
34112PRO	230 cm	7 Fr	12 mm	13.5 mm	15 mm	55 mm
34115PRO	230 cm	7 Fr	15 mm	16.5 mm	18 mm	55 mm
34118PRO	230 cm	7 Fr	18 mm	19 mm	20 mm	55 mm
Included in delivery per package: 1						

Biliary Balloon Dilator

Part number	Catheter length	Catheter Ø	Balloon Ø			Balloon length
3410630PRO	200 cm	7 Fr	6 mm	7 mm	8 mm	30 mm
3410830PRO	200 cm	7 Fr	8 mm	9 mm	10 mm	30 mm
3411030PRO	200 cm	7 Fr	10 mm	11 mm	12 mm	30 mm
3411230PRO	200 cm	7 Fr	12 mm	13.5 mm	15 mm	30 mm
Included in delivery per package: 1						