

Polypectomy Loops

Twistable polypectomy loops for single use

- 2.3 mm \varnothing incl. handle
- For working channels **2.8 mm in diameter and larger**
- **The following designs are available:**
 - Hexagonal, multifilament
 - Asymmetric, multifilament
 - Oval, multifilament
 - Oval, monofilament



Part number	Form	Loop opening	Length	Loop wire
NOE372216-G	Hexagonal	25 mm	230 cm	multifilament
NOE372217-G	Hexagonal	35 mm	230 cm	multifilament
NOE3222-G	Asymmetric	15 mm	230 cm	multifilament
NOE3322-G	Asymmetric	25 mm	230 cm	multifilament
NOE342214-G	Oval	15 mm	230 cm	multifilament
NOE342216-G	Oval	25 mm	230 cm	multifilament
NOE342217-G	Oval	35 mm	230 cm	multifilament
NOE342212M-G	Oval	6 mm	230 cm	monofilament
NOE342213M-G	Oval	10 mm	230 cm	monofilament
NOE342214M-G	Oval	15 mm	230 cm	monofilament
NOE342216M-G	Oval	25 mm	230 cm	monofilament
NOE342217M-G	Oval	35 mm	230 cm	monofilament

Included in delivery per package: 5 pieces