

Sclerotherapy Needle

For single use

The Sclerotherapy Needle has a metal head that ensures controlled advancement. The needle can be used for all working channels with a diameter of 2.8 mm and larger.



- **2.3 mm Ø**
- For working channels **2.8 mm in diameter and larger**
- **With a metal head** for controlled advancement of the needle

Part number	NET2522-C5	NET2522-G5
Needle length	5 mm	5 mm
Needle diameter	0.7 mm	0.7 mm
Total length	180 cm	230 cm
Included in delivery per package: 5 pieces		