

Duo System Premium

For convenient single-hand operation

The DUO System PREMIUM is a correctable localisation system for mammography. It offers doctors and patients top quality, absolute reliability and the highest degree of comfort in the preoperative marking of suspicious tissue.

The marker wire of the DUO System PREMIUM can be easily retracted into the cannula to correct its position thanks to its special form and then be released again.

The double arches provide firm anchorage in the tissue and thereby prevent dislocations in any direction. The depth position of the marked lesion can be detected based on the various properties of the wire.

DUO System PREMIUM offers easy and precise handling as well as a gentle application and ensures optimum visibility in ultrasound.

- **Single-handed operation**
- **Correctable marker wire**
- **By means of DUO arches** firm anchorage in the tissue
- **Sharp wire tips** for hard lesions
- **Highly flexible** marker wire
- **Depth orientation** at the cannula and wire
- **Ultrasound marking and ultrasound cut** at the cannula tip
- **Fixation piece** allows for the secure preoperative fixation of the wire



Part number	271770	271771	271773
Gauge	20 G	20 G	20 G
Diameter	0.95 mm	0.95 mm	0.95 mm
Length	50 mm	90 mm	120 mm
Included in delivery per package: 10 pieces			