

Duo System

Correctable localisation system

The Duo System is a correctable localisation set for preoperative marking of non-palpable, suspicious mammary lesions. It makes the intra-operative localisation of findings easier for the surgeon. For required positional correction, the marker wire can easily be retracted back into the cannula and released again for marking.

- **Correctable marker wire**
- **Double-loop technique** for secure anchoring
- **Highly flexible** wire
- **Atraumatic, extremely sharp** cannula cut
- **Fixation piece** to prevent distal dislocation



Duo System

Part number	271440	271441	271443
Gauge	20 G	20 G	20 G
Diameter	0.95 mm	0.95 mm	0.95 mm
Length	50 mm	90 mm	120 mm
Included in delivery per package: 10 pieces			

MR Duo System

Part number	601609	601611
Gauge	20 G	20 G
Diameter	0.95 mm	0.95 mm
Length	90 mm	120 mm
Included in delivery per package: 10 pieces		