


Tumark clip marker

Gentle minimally invasive marking of tumours in the breast

The marker of the Tumark line can pinpoint markings in the soft tissue through the placement of a marker before or during neoadjuvant chemotherapy. Furthermore, the system is ideal for the marking of biopsy collection sites. The ergonomic handle of the clip-marking system enables single-handed operation to the user during the simultaneous guiding of the ultrasonic head during the intervention.



- **Precise, safe and gentle** clip marking of soft tissue
- **Excellent visibility** in ultrasound imaging
- **Firmly fixed anchorage** in the tissue

Part number	Description	Cannula		
		Diameter	Gauge	Length
271500	Tumark Q	1,20 mm	18 G	120 mm
271551	Tumark Professional X-Shape	1,20 mm	18 G	120 mm
271560	Tumark Professional	1,20 mm	18 G	120 mm
271589	Tumark Vision	1,20 mm	18 G	100 mm
271590	Tumark Vision	1,20 mm	18 G	120 mm
601570	Tumark MRI	1,20 mm	18 G	120 mm
601590	Tumark Vision MRI	1,20 mm	18 G	120 mm

Included in delivery per package: 10 pieces

For other sizes, please contact our product specialist