


VenoSeal™

Relief for venous leak

For a man, venous leak means that the veins in the penis cannot hold the blood for long enough during an erection, resulting in venous reflux, usually during sexual intercourse (activity). This means that getting an erection is not a problem, but maintaining it is. VenoSeal™ supports you in this.



All benefits at a glance:

- Easy to use
- Comfortable to wear thanks to the soft material
- Very discreet thanks to its colour
- Can be used in addition to PDE5 medications
- Reusable (not disposable)

Brief instructions at a glance:

- Position the loop at the base of the penis (before or during an erection)
- Make sure that the end with the ball (**red arrow**) is facing upwards
- Hold the end without the ball (**blue arrow**) between the thumb and the index finger
- Pull in a horizontal direction away from the ball until the ring is firmly secured
- Only pull tight enough to ensure that blood reflux is prevented (minimal pressure)
- After 30 minutes, the ring should be removed by pulling on the end with the ball and thereby widening the loop

Important information:

- VenoSeal™ is not intended to replace a visit to your doctor!
- Erectile dysfunction could be an indication of a variety of conditions